

PERCEPTRON-XFMS: An Open-Source Web Server for Analysis of X-ray Footprinting with Mass Spectrometry Data

User Guidelines

**Biomedical Informatics and Engineering Research Laboratory,
Department of Life Sciences,
Lahore University of Management Sciences,
Pakistan.**

biolabs.lums.edu.pk/BIRL

Last Modified: 19-03-2025


1		
2	Contents	
3	1. Step-by-Step Guidelines for using MASCOT	3
4	A. Mgf file for MASCOT	3
5	B. Generating Mascot Result File	4
6	2. Preparation of Input Files	5
7	3. Step-by-Step Guidelines for Using Webserver.....	7
8	A. Downloadable Result Files	10
9		
10		

1. Step-by-Step Guidelines for using MASCOT

PERCEPTRON-XFMS is a highly interactive web server developed to assist in XFMS data analysis. The pipeline initiates with input of LC-MS/MS data along with parameter files. Upon receiving the dose-dependent LC-MS/MS data, the user must perform peptides identification process using MASCOT. The step-by-step process for CheY protein is listed below.

A. Mgf file for MASCOT

1. Open MSConvert GUI and browse the MS/MS file with the highest dosage.
2. Select the output format as '.mgf' and press the 'Start' button.



B. Generating Mascot Result File

1. Open Mascot MS/MS Ions search through the following link and select: UP625_E_Coli_K12 (AA) database
2. Select the desired variable modification.
3. Upload the '.mgf' file.
4. Start searching and download the results.

MASCOT MS/MS Ions Search

Your name **Email**

Search title

Database(s) UP625_E_coli_K12 (AA) 1

Taxonomy

Enzyme **Allow up to** missed cleavages

Quantitation

Crosslinking

Fixed modifications

Variable modifications 2

☒ Display all modifications

Peptide tol. \pm Da **# ^{13}C**

Peptide charge

Data file Choose File STB4edited.mgf 3

Data format

Instrument

☐ Decoy

MS/MS tol. \pm Da **Monoisotopic** ☒ **Average** ☐

Precursor m/z

☐ Error tolerant

Target PSM FDR

Start Search ... 4

2. Preparation of Input Files

Create 2 separate folders for “MS/MS Input Files” containing (1) experimental data, and (2) “Additional Input Files” having the following parameters.

I. **Folder 1 (MS/MS Input Files)** - This folder contains the LC-MS/MS mzXML experimental files.

II. **Folder 2 (Additional Input Files)** - This folder contains the following files.


a) **Mascot.xlsx** – Mascot result file

b) **Protein_structure.pdb** – protein structure file

c) **SASA.xlsx** – Solvent accessible surface area (SASA) file generated by POPS server. The file includes detailed residue SASA information including Residue Name, Chain, Residue Number, Hydrophobic SASA (Phob/A^2), Hydrophilic SASA (Phil/A^2), Total SASA (Total/A^2), Relative SASA ($\text{Q}(\text{SASA})$), Number of Overlaps ($\text{N}(\text{overlap})$), and Surface Area (Surf/A^2). Currently, PERCEPTRON-XFMS plots the $\log(\text{PF})$ value against the Total SASA by default.


d) **DoseResponseInformation.txt** – A text file that contains the X-ray doses used in the study.

PERCEPTRON-XFMS: User Guidelines




1

- 1 **Step 3** - Upload the LC-MS/MS experimental data. PERCEPTRON-XFMS requires experimental
- 2 mzXML-formatted data to be arranged as three separate replicates for each X-ray dose.
- 3 Specifically, the user must provide three replicates per dose, with each replicate stored as an
- 4 individual mzXML file. The total number of mzXML files required is $n \times 3$, where n is the number
- 5 of doses used in the experiment. Names of files must strictly follow the same naming structure as
- 6 the DoseResponseInfo file i.e., “Dose0.mzXML” etc.



- 7
- 8 **Step 4** – Protein Contact Analyses (Advanced Users)
- 9 To analyze protein-water hydrogen bond graphs click the ‘Bridge Algorithm’ box.
- 10 Click on the ‘Frustratometer Algorithm’ box for additional analysis of protein residue energetics.

PERCEPTRON-XFMS: User Guidelines




1

2 Submit the Query Search by clicking on the ‘Submit’ button.

3 **Step 5-** Click on ‘Search Results & History’ to see the search results. Click on your “Search


4 Title” row.



5

6 **Step 6** – Download results by clicking ‘**Results Download**’ button option provided on top right.

7 The results will be downloaded as a *.zip file.



1

2 A. Downloadable Result Files

3 Downloaded result (*.zip) file contains multiple following subfolders and files.

4 i. Several subfolders correspond to different peptides in the protein and each folder contains:

5 a) Subfolders containing rate of oxidation for individual replicate.

6 b) Dose-response plots of the rate of oxidation.

7 c) **Mean_R1_R2_R3.xls** file contains the mean values of (1 - Rate of oxidation) for each
8 modified residue reported in that peptide.

9 ii. **PeptideInfo.xls** file which includes information regarding the peptide sequence and their
10 corresponding peptide folder.

11 iii. **PF_SASA_tab.xls** - contains a table with log protection factor of each modified residue.

12 iv. **SASAmain.png** - a graphical representation of the protection factor data with SASA values on
13 the x-axis and log (PF) on the y-axis.

14 v. **Modified.pdb** - a pdb file incorporated with log (PF) values.

1 vi. **ResultsBridge.xlsx** - an Excel file that contains centrality scores of the residues using the Bridge
2 algorithm.

3 vii. **FrustratometeR-Results** include images that quantify the degree of local frustration in each
4 protein. Users also get tables with the calculated Frustration indexes and scripts to locally explore
5 and manipulate the structures.

6

7